CARDIO

EC038M

CE 08/12/2016

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The compression garment facilitates wound healing after surgery. It is intended to support the tissues by maintaining them in their new shape, thereby reducing post-operative edema and limiting skin laxity. Thermoregulating Coolmax® fabric allows moisture to be evacuated towards the outer layer of the garment.

Application notes

For optimal effectiveness, the compression garment should be applied immediately after surgery. The duration of the treatment depends on the medical prescription. The first application should be overseen by a health professional. The compression garment must not be worn directly on injured skin.

The sizes indicated correspond to the measurements taken before surgery. Other sizes can be made to measure.

Care

Wash your garment at 30° C maximum with a natural soap-based detergent. Do not use a tumble dryer. We recommend to wash the garment before wearing it for the first time.

Side effects

In very rare cases, skin intolerances or allergies may occur. In this case, consult your physician. In the event of a serious incident with this medical product, please report it to the manufacturer and to the competent authorities of the country you live in.

Contraindications

Don't use the compression garment if you detect an allergy or skin intolerance. Seek advice from your physician!